


MyoStim®

TEMPORARY MYOCARD-ELECTRODES



Temporary Myocard-Electrodes


MyoStim®

In combination with an external pacemaker, myocardial electrodes serve to temporarily stimulate the heart during and after cardiac surgery. Several variants are available:

Number of Poles	Unipolar, Bipolar (bifurcated/ in-line) or Quadripolar (bifurcated/ in-line)
Fixation on the Heart	Tines, Zigzag, V-Hook, Loop Tines+ V-Hook, Tines + Loop <small>(only Quadripolar)</small>
Length of the Myocard-Electrode	60 cm or 220 cm
Pacemaker Adapter	Myoconnector, Pin or Breaking Needle
Color Code	White, Blue or White/Blue <small>(only Quadripolar)</small>

Other variants on request

MyoStim Unipolar




Product Code for Unipolar:

HW.	5	0		C			
Length:				Color:	Fixation:	Adapter:	
060 = 60 cm				empty = Blue	TIN = Tines	empty = Breaking Needle	
220 = 220 cm				W = White	TZ = Zigzag	M = Myoconnector	
					TV = V-Hook		
					TL = Loop		

Example: HW.50060CTINM – Unipolar, 60 cm, Blue, Tines, Myoconnector
 HW.50220CWTV – Unipolar, 220 cm, White, V-Hook, Breaking Needle

Packaging unit: 25 pieces

MyoStim Bipolar bifurcated



Product Code for Bipolar bifurcated:


HW.	5	2		C			
Length:				Color:	Fixation:	Adapter:	
060 = 60 cm				W = White	TIN = Tines	empty = Pin	
220 = 220 cm				B = Blue	TZ = Zigzag	M = Myoconnector	
					TV = V-Hook		
					TL = Loop		

Example: HW.52060CWTZ – Bipolar bifurcated, 60 cm, White, Zigzag, Pin

HW.52220CBTVM – Bipolar bifurcated, 220 cm, Blue, V-Hook, Myoconnector

Packaging unit: 25 pieces

MyoStim Bipolar in-line



Product Code for Bipolar in-line:

HW.	5	3		C			
Length:				Color:	Fixation:	Adapter:	
060 = 60 cm				W = White	TIN = Tines	empty = Pin	
220 = 220 cm				B = Blue	TZ = Zigzag	M = Myoconnector	
					TV = V-Hook		
					TL = Loop		


Example: HW.53060CBTINM – Bipolar in-line, 60 cm, Blue, Tines, Myoconnector

HW.53220CWTZ – Bipolar in-line, 220 cm, White, Zigzag, Pin

Packaging unit: 25 pieces

Temporary Myocard-Electrodes

MyoStim Quadripolar bifurcated



Product Code for Quadripolar bifurcated:


HW.	5	4		C		
Length:			Fixation:		Adapter:	
060 = 60 cm			TIN = Tines		empty = Pin	
220 = 220 cm			TZ = Zigzag		M = Myoconnector	
			TV = V-Hook			
			VTIN = Tines + V-Hook			
			TINL = Tines + Loop			

Example: HW.54060CTIN – Quadripolar bifurcated, 60 cm, With/Blue, Tines, Pin

Packaging unit: 20 pieces

HW.54220CTVM – Quadripolar bifurcated, 220 cm, White/Blue, V-Hook, Myoconnector

MyoStim Quadripolar in-line



Product Code for Quadripolar in-line:

HW.	5	5		C		
Length:			Fixation:		Adapter:	
060 = 60 cm			TIN = Tines		empty = Pin	
220 = 220 cm			TZ = Zigzag		M = Myoconnector	
			TV = V-Hook			
			VTIN = Tines + V-Hook			
			TINL = Tines + Loop			

Example: HW.55060CTIN – Quadripolar in-line, 60 cm, White/Blue, Tines, Pin

Packaging unit: 20 pieces

HW.55220CTVM – Quadripolar in-line, 220 cm, White/Blue, V-Hook, Myoconnector

Osypka Medical GmbH

Albert-Einstein-Strasse 3
12489 Berlin, Germany
Phone: +49 (30) 6392 8300
Fax: +49 (30) 6392 8301
E-Mail: mail@osypkamed.com
www.osypkamed.com



MyoStim is a trademark of Osypka Medical
Copyright © 2019 Osypka Medical

