

• • • •

Trolley and Medical Grade Tablet PC


In this setup, you may install a printer on the printer shelf, and the customer may print reports to include them in the patient file.





Specifcations:

- 1- Hemoview Screen (8 bars + Signal Waveforms).
- 2- Trend Screen with labels for custom made time markers and export options.
- 3- Report Screen with creating Comparative Time Reports.
- 4- Passive Leg Raising Test screen (To test patient response for fluids-very important test in ICU-).
- 5- Optional Management Screen.
- 6- Unlimited reports could be created for each patient.
- 7- Unlimited patient data could be stored and kept on the iControl Software.
- 8-Patient data storage could include waveforms.
- 9- Powerful portable workstation.
- 10- Excellent printing and storage capabilities.
- 11- Data export to Hospital Information Systems HIS.
- 12- Perfect for intra-hospital use.



Discription: Trolley and Medical Grade tablet PC for ICON and ICONCORE





